

**National Initiative on Undergraduate Science - Teacher Development
Workshop
(NIUS-TD Workshop)
November 3-8, 2019**

**Homi Bhabha Centre for Science Education
Tata Institute of Fundamental Research**

In collaboration with

Gogate Joglekar College, Ratnagiri

Organic Chemistry Task IB






This workshop is organized by Centre for Excellence in Science and Mathematics Education (CESME-HBCSE) and set up under the Pandit Madan Mohan Malaviya National Mission for Teachers and Teaching Programme

This experiment has been adapted from 45th International Chemistry Olympiad held in Moscow, Russia in 2013 which is available in the public domain and is accessed from: <http://www.icho2013.chem.msu.ru/en/problems-solutions/icho-2013-problems>

Tracking the supplied samples of substituted aldehydes by hydrazone derivativesList of Glassware and equipment

| Item | Quantity |
|---|----------|
| 5 mL Plastic tube with screw neck labeled "1" | 1 |
| 5 mL Plastic tube with screw neck labeled "2" | 1 |
| Burette stand | 1 |
| 50 mL beaker | 2 |
| 15 mL test tube | 4 |
| Spatula | 2 |
| 500 mL plastic washer bottle | 1 |
| 10 mL measuring cylinder | 1 |
| Glass rod | 1 |

List of Chemicals

| Reagent | Quantity | Labeled | Safety |
|--|----------------------|---|---|
| 2,4-Dinitrophenylhydrazine | 200 mg each, 2 vials | 2,4 DNP |  |
| Sulfuric acid, concentrated | 1 mL each, 2 tubes | H ₂ SO ₄ concentrated |  |
| Aldehyde solution 1 mmol in ethanol | 4 mL each, 2 bottles | Aldehyde 1 and Aldehyde 2 |  |
| NaOH solution | 10 mL | NaOH 2M |  |
| Acetone | 10 mL | Acetone |  |
| NaHCO ₃ solution | 10 mL | NaHCO ₃ 10% | |
| Ethanol | | | |

On the table for the common use

Ethanol, Liquid waste container, Solid waste container, Distilled Water, Filter paper, pH papers

Hazard Symbols

| | | | | | |
|--------------|---|-----------|---|-----------------------|---|
| Corrosive |  | Flammable |  | Irritant |  |
| Carcinogenic |  | Toxic |  | Hazard to environment |  |


Tracking the supplied samples of substituted aldehydes by hydrazone derivatives**Introduction:**

Hydrazone are organic compounds with the general formula, $R_1R_2C=N-NHR_3$. Reaction of hydrazine or substituted hydrazines (like 2,4-dinitrophenylhydrazine) with ketones or aldehydes lead to formation of hydrazones, (imine derivatives). The properties of hydrazones are largely due to the presence of the $>C=N-N<$ structural unit, which contains two nitrogen atoms that are nucleophilic in nature. The carbon atom may act as either electrophile or nucleophile according to the reaction conditions and environment. The $C=N$ bond of dialkyl hydrazones can be hydrolyzed, oxidized and reduced that gives variety of products. Similarly the $N-N$ bond can be reduced to the free amine.

These derivatives are highly coloured solids, often stable and crystalline. You must have used formation of hydrazones as a qualitative test for confirming the presence of a carbonyl group. In fact, the concentration of low molecular weight aldehydes and ketones can be determined by converting them to their hydrazone derivatives. Hydrazone derivatives have various important pharmacological activities and are used for synthesis of a wide variety of medicinally active compounds.

References: i) Reactions, Mechanisms, and Structure, John Wiley & Sons, Milton, Australia, 6th edition, 2007.
ii) <http://dx.doi.org/10.1155/2014/761030>

In the current task, you will be preparing substituted hydrazones using 2,4-dinitrophenylhydrazine for the given samples of two substituted benzaldehydes (shown below). In fact, using the colour and some simple qualitative tests, you will identify the substituted benzaldehydes given to you.



Procedure**Preparation of 2,4-dinitrophenylhydrazones**

Weigh 0.200 g of given 2,4-dinitrophenylhydrazine and transfer it to a 50 mL beaker and place a glass rod inside the beaker. Add 2 mL of the given concentrated sulfuric acid to this beaker in a dropwise manner using a **dropper and with constant stirring**. Then, using measuring cylinders, add 1.6 mL of water and 4 mL of ethanol to the reaction mixture. At the end, add one of the aldehyde sample solution and note its label (either **sample 1** or **sample 2**). A bright coloured precipitate will start forming at once. Continue stirring with the glass rod for 10 min, then add 10 mL of water and stir for another 3 min.

Separation and purification of the product

Filter the product using the Büchner funnel under suction. Wash the precipitate with small amounts of water (2 -3 mL each time) until the pH of the filtrate is almost neutral. Lastly, wash the precipitate twice with ethanol (2-3 mL each time). Air-dry the solid on the filter paper and record the mass of the product in the answer sheet.

Repeat the entire procedure with the other sample of the aldehyde provided to you.

Answer Sheet

| | | | |
|---|--|--|--|
| Plastic Vial/stoppered tube 1 (Sample 1) | | Plastic Vial/ stoppered tube 2 (Sample 2) | |
| Mass of empty vial/ stoppered tube | | Mass of empty vial/ stoppered tube | |
| Mass of vial + product | | Mass of vial + product | |
| Mass of product | | Mass of product | |

1. The structures of the aldehydes are given to you. Draw the structures of the two hydrazones that will be formed.

| Product 1 (sample 1) | Product 2 (sample 2) |
|-------------------------|-------------------------|
| | |

- 2a) In your opinion, will the hydrazones formed exhibit stereoisomerism?

Hydrazone formed for Sample 1 Yes No

Hydrazone formed for Sample 2 Yes No

- b) If your answer is yes for one or both, the kind of stereoisomerism will be

3. In your opinion, whether the acidic pH of the reaction mixture is a crucial factor which will affect the yield of the above synthesis? Explain your answer in brief.

4. In your opinion, whether the synthesis reaction will proceed in

- | | | | | |
|------------|-----|--------------------------|----|--------------------------|
| a) neutral | Yes | <input type="checkbox"/> | No | <input type="checkbox"/> |
| b) basic | Yes | <input type="checkbox"/> | No | <input type="checkbox"/> |

5. If the above synthesis is carried out with **p-nitro benzaldehyde** in acidic medium the rate of formation of product will be:

- | | |
|--|--------------------------|
| a) slower as compared to the synthesis done by you | <input type="checkbox"/> |
| b) faster | <input type="checkbox"/> |
| c) no change | <input type="checkbox"/> |
| d) cannot say as more information will be required | <input type="checkbox"/> |

Characterization

Dissolve a **single crystal** of each hydrazone in a test tube in 10 mL of acetone and divide into two parts.

6. a) To I part add 1 drop of **saturated** NaHCO_3 solution and shake the resultant solution well. Record your observations of the solutions color changes in the box.

Observation for 2,4-DNPH derivative of Sample 1

The color of the solution does not change

The color changes significantly

Observation for 2,4-DNPH derivative of Sample 2

The color of the solution does not change

The color changes significantly

b) To the second part of the solution, add **2 drops** of NaOH solution and shake the resultant solution well. Record your observations of the solutions color change in the box.

Observation for 2,4-DNPH derivative of Sample 1

The color of the solution does not change

The color changes significantly

Observation for 2,4-DNPH derivative of Sample 2

The color of the solution does not change

The color changes significantly

7. What structural features of your products explain the color change (if any) in the reaction with NaHCO_3 ? Tick the appropriate box/es.

- presence of MeO group at position 4 in the benzene ring;

presence of the OH group at position 2 in the benzene ring;

presence of both $>\text{C}=\text{N}-\text{N}<$ unit

presence of both $>\text{C}=\text{N}-\text{N}<$ and OH

8. What structural features of your products explain the color change (if any) in the reaction with NaOH ? Tick the appropriate box/es.

- presence of MeO group at position 4 in the benzene ring;

presence of the OH group at position 2 in the benzene ring;

presence of both $>\text{C}=\text{N}-\text{N}<$ unit

presence of both $>\text{C}=\text{N}-\text{N}<$ and OH

9. Which of the listed processes is responsible for the color change (if any) observed in the reaction of 2,4-dinitrophenylhydrazones with aqueous NaOH ? Tick the appropriate box.

- alkaline hydrolysis dehydration hydration

deprotonation dehydrogenation

10. Draw the structures of the main organic species present in each test reaction medium in the answer box below.

| Initial aldehyde | Solution with NaHCO_3 | Solution with NaOH |
|------------------|--------------------------------|-----------------------------|
| Sample 1 | | |
| Sample 2 | | |

11. Calculations

| | Mass of aldehyde | Molar Mass of aldehyde | Molar Mass of hydrazone | Mass of hydrazone |
|----------|------------------|------------------------|-------------------------|-------------------|
| Sample 1 | | | | |
| Sample 2 | | | | |

12. Calculate the theoretical yield and the percentage yields of the hydrazones.

| |
|----------|
| Sample 1 |
| Sample 2 |